

M9111 Filtration Assemblies

Filtration Assemblies are designed for high pressure liquid chromatography (HPLC), particulate, and microbiological contamination filtrations.

Sintered glass support provides even, flat support structure for membrane filters. They are available in 47-mm or 90-mm filter sizes. Filtering flasks are also available, but sold separately.



Features

- *Coarse fritted glass support base*
- *Graduated glass funnel with 300ml or 1,000ml capacity*
- *Anodized aluminum clamp*
- *#8 silicone stopper with 9/16-inch hole*

Specifications

Size:	90 mm
Diameter:	90 mm
Filter Sizes:	47mm/ 300ml capacity 90mm/ 1,000ml capacity